



DO LESS.

CARE MORE.



Clear Direction.  From Diagnosis to Care.

**medrad**® Centargo  
CT Injection System

# DO LESS.

**Centargo** minimizes your injector touch time by automating manual processes.

# CARE MORE.

**Centargo** helps you focus on what matters most in every CT scan – the patient.



# With Efficiency.

Centargo is designed for a simple, and automated workflow, easy daily set-up and quick change time between scans.\*

Centargo's daily set-up process can be completed in under **2 minutes**.



The simple **snap-in** patient line **auto-primes** upon insertion and is ready for the next patient in less than **20 seconds**.



# With Confidence.

Automated data entry and the intuitive user interface help you confidently complete your daily tasks.

**Touch screen** within the scan room so that you can be close to your patient.



The **bar code reader** reduces manual data entry and provides easy traceability and access to contrast and injection details.\*\*



\* Centargo Operation Manual, Centargo System Setup Timing DF-1502305

\*\* When combined with Automated Documentation software

# With Centargo.

## For easy set-up and patient-to-patient change over:

- › Pre-assembled 24-hour Day Set with optional replaceable spikes
- › Fast changeover using snap-in Patient Line

## Designed with safety in mind:

- › Inlet and outlet air detection in addition to automatic air removal within reservoirs
- › Recessed connection on the Day Set
- › Two check valves in the Patient Line

## Together with performance:

- › Piston-based fluid delivery
- › Heat maintainers for contrast reservoirs

## And convenience:

- › Integrated prime container
- › Dual screen configuration within scan and control room
- › Training videos on Control Room Unit
- › VirtualCare for remote monitoring, updates, and diagnostics

# With Workflow Solutions.

- › Integrated barcode reader for contrast data capture
- › Automated Documentation for access to contrast and injection information in PACS
- › Smart Protocols to calculate individualized injection protocols
- › Advanced Scanner Connectivity for protocol exchange and device synchronization



## With Quality.

### Centargo is built on the trusted quality of Bayer MEDRAD®.

- › Piston-based fluid delivery and personalised contrast media protocol software.\*
- › Three inlet air sensors and one outlet air sensor as well as reservoir air detection and automatic air removal.
- › Cordless mobility through a battery powered, Wi-Fi enabled configuration.

\*P3T® and Protocol Assistant Software

Think ahead. Plan ahead. Stay ahead.  
Bayer is your partner for peace  
of mind and patient care.

## With Bayer.

Bayer relies on our decades of experience, comprehensive service networks, expert clinical support and dedicated sales teams to make sure that you receive the solutions you need to improve your efficiency, confidence and quality.



# System Specifications\*

Description	Specification
Volume Range	1 mL to 200 mL for contrast and flush phases; 1 mL to 400 mL for DualFlow phases (in 1 mL increments)
Flow Rate Range	0.1 to 10 mL/sec in 0.1 mL/sec increments
DualFlow	Simultaneous injection of contrast and saline in increments of 5%
Patency Check / Test Inject	Real time flow rate adjustment and multi-phase test injections
Programmable Pressure Limit	From 50 psi to 300 psi in 1 psi increments (also displayable in kPa and kg/cm <sup>2</sup> )
Programmable Delay	1 second to 900 seconds in 1 second increments
Pause	Maximum pause time of 20 minutes
Fluid Loading	Contrast: 50 mL – 500 mL Saline: 50 mL to 1000 mL
Fluid Reservoirs (volume)	Three 200 mL reservoirs (two contrast and one saline)
Injection Capabilities	Maximum of 6 phases per injection. Maximum of 10 injections per exam
Fluid Delivery Method	Piston-based
Day Set	Pre-assembled, labeled for 24-hour use, designed for set-up in less than 2 minutes
Patient Line	250 cm, single hand design, two check valves, auto-primers upon insertion
Displays	Two independent, Wi-Fi enabled, full color, touch screen units with modern graphical user interfaces located in the scan room and control room
Air Detection	Three inlet air sensors and one outlet air sensor Includes reservoir air detection and automatic air removal
Battery Performance	Up to 16 hours or up to 60 injections
Integrated Prime Container (vol)	680 mL
Injector Dimensions (WxDxH)	76.2 x 68.6 x 167.6 cm
Injector Weight (approx.)	76 kg (Pedestal configuration)
Available Scanner Interface	ISI 700, ISI 900 and Connect.CT
Smart Protocols	P3T <sup>®</sup> algorithms calculate individualized injection protocols for cardiac CT, CT pulmonary angiography and CT imaging of abdominal organs Protocol Assistant provides basic protocol calculations for a range of procedure types
eGFR Calculator	Available on Control Room and Scan Room Units

\* Centargo Operation Manual, Centargo System Setup Timing DF-150205

## Ordering Information

### Injector Configurations

➤ Centargo Pedestal System with Batteries **CENT-SYS-BAT**

### Disposables

➤ Centargo Day Set (24-hour use) **CENT-DS**  
 ➤ Centargo Patient Line **CENT-PL**  
 ➤ Centargo Replacement Spikes (optional) **CENT-RS**

### Smart Protocols

➤ P3T Personalized Protocols **CENT-P3T-PAT**  
 ➤ P3T Protocol Assistant **CENT-PAT**  
 ➤ PAT Protocol Assistant **CENT-P3T**

### Automated Documentation

➤ Modality Worklist **CENT-AUTO-DOC**  
 ➤ PACS Interface **CENT-MWL**  
**CENT-PACS**

### Processing Productivity<sup>§</sup>

➤ RIS Interface **CENT-RIS**  
 ➤ Dictation Interface **CENT-SR**  
 ➤ Radimetrics Interface **CENT-REP**

### Scanner Connectivity<sup>†</sup>

➤ ISI 900 / ISI 700 **ISI2**  
 ➤ Connect.CT **CENT-ISI-CCT**

§ With Additional Purchase

† Only available via OEM



2920 Matheson Blvd East  
 Mississauga ON L4W 5R6  
 Phone: (800) 268-1432  
 Fax: (800) 567-1710

Bayer reserves the right to modify the specifications and features described herein or to discontinue any product or service identified in this publication at any time without prior notice or obligation. Please contact your authorised Bayer representative for the most current information.

All patient data that appear in this document are fictitious. No actual patient information is shown.

Bayer Cross, MEDRAD, MEDRAD Centargo, Centargo, P3T, VirtualCare and Radimetrics are trademarks owned by and/or registered to Bayer in the U.S. and/or other countries. Other trademarks and company names mentioned herein are properties of their respective owners and are used herein solely for informational purposes. No relationship or endorsement should be inferred or implied.

CEN001E ©2025 Bayer. This material may not be reproduced, displayed, modified or distributed without the express prior written consent of Bayer.